

## CONTACT LENS PARAMETER CHEAT SHEET

<b>Manufacturer/Brand</b>	<b>BC</b>	<b>DIA</b>
<b><u>Johnson &amp; Johnson-(Vistakon) (2 WK lenses)</u></b>		
Acuvue 2	8.3, 8.7	14.0
Acuvue 2 Colours	8.3	14.0
Acuvue Advance Plus	8.3, 8.7	14.0
Acuvue Advance for Astigmatism	8.6	14.5
Acuvue Bifocal	8.5	14.2
Acuvue Oasys	8.4, 8.8	14.0
Acuvue Oasys for Astigmatism	8.6	14.5
Acuvue Oasys for Presbyopia	8.4	14.3
Acuvue One Day Moist(day)	8.5, 9.0	14.2
<b><u>Cooper Vision</u></b>		
Avaira (2 wk)	8.4, 8.5	14.2
Biofinity (mo)	8.6	14.0
Biofinity Multifocal (mo)	8.6	14.0
Biofinity Toric (mo)	8.7	14.5
Frequency 55 (mo)	9.0, 8.7, 8.4	14.2
Frequency 55 Aspheric (mo)	8.7, 8.4	14.4
Frequency 55 Multifocal(mo)	8.7	14.4
Frequency 55 Toric & Toric XR (mo)	8.7, 8.4	14.4
Proclear (Biomedics) EP (2 WK)	8.7	14.4
Proclear Multifocal (mo)	8.7	14.4
Proclear (mo)	8.6, 8.2	14.2
Proclear Toric (mo)	8.4	14.4
Proclear Toric XR (mo)	8.7, 8.4	14.4
<b><u>Alcon-(Ciba Vision)</u></b>		
Air Optix Aqua (mo)	8.6	14.2
Air Optix Aqua Multifocal (mo)	8.6	14.2
Air Optix for Astigmatism (mo)	8.7	14.5
Air Optix Night & Day (mo)	8.4, 8.6	13.8
Freshlook Visitint & Colors (2wk)	8.6	14.5
Freshlook Toric (2wk)	8.4	14.5
O2 Optix (mo)	8.6	14.2
<b><u>Bausch &amp; Lomb</u></b>		
Purevision (mo)	8.6, 8.3	14.0
Purevision 2 HD(mo)	8.6	14.0
Purevision Toric (mo)	8.7	14.0
Purevision Multifocal (mo)	8.6	14.0
Soflens 38 (2wk)	9.0, 8.7, 8.4	14.0
Soflens Multifocal (2wk)	8.5, 8.8	14.5
Soflens Daily Disposable (day)	8.6	14.2
Soflens Toric (2wk)	8.5	14.5